



## CHATA Biosystems Caffeine Linearity Standards

Chata Biosystems offers NIST traceable Caffeine Linearity Standards pre-packaged in HPLC autosampler vials. Just pop them in your autosampler and you are ready to go!

Standards are provided in duplicate sets. Each set contains a water blank and vials, 1.8 ml each, in 0.5, 1.0, 5.0, 25.0 and 50.0 mg/L concentrations, a total of twelve vials per package. Supporting stability data demonstrate a 6-month expiry.

- **Increased Safety**  
-no glass ampoules to break open
- **No Liquid Transfer**  
-the vials are autosampler ready
- **No Syringe Cleaning**  
-extraneous glassware is unneeded
- **Autosampler Ready**  
-unsurpassed ease of use
- **Significant Time Savings**

These standards are typically used for IQ/OQ and validation testing of HPLC systems. They are fully traceable and are manufactured under cGMP.

[http://208.101.44.196/uploads/ChataCHEMNECT\\_Dissoluton\\_Media.pdf](http://208.101.44.196/uploads/ChataCHEMNECT_Dissoluton_Media.pdf)