



CHATA Biosystems Caffeine Linearity Standards

Chata Biosystems offers NIST traceable Caffeine Linearity Standards pre-packaged in HPLC autosampler vials. Just pop them in your autosampler and you are ready to go!

Standards are provided in duplicate sets. Each set contains a water blank and vials, 1.8 ml each, in 0.5, 1.0, 5.0, 25.0 and 50.0 mg/L concentrations, a total of twelve vials per package. Supporting stability data demonstrate a 6-month expiry.

- **Increased Safety**
 - no glass ampules to break open
- **No Liquid Transfer**
 - the vials are autosampler ready
- **No Syringe Cleaning**
 - extraneous glassware is unneeded
- **Autosampler Ready**
 - unsurpassed ease of use
- **Significant Time Savings**



These standards are typically used for IQ/OQ and validation testing of HPLC systems. They are fully traceable and are manufactured under cGMP.

Reference Part # 15015

Chata Biosystems, Inc

323 S College Ave
#5 Clocktower Square
Fort Collins, CO 80524
Telephone 877-248-2428 |
Fax 970-221-9488
Email: sales@chatasolutions.com
www.chatasolutions.com